

# Common Reporting Tool 1.3 (CRT)

GENERATE YOUR AVS AND VALPROBE RT  
REPORTS FROM YOUR PC



The Common Reporting Tool software is a comprehensive reporting utility that can be installed on a separate Windows PC and allows generating reports from Kaye Validator AVS or ValProbe RT Qualification study files to document the results of your validation studies. Qualification study files, which have been generated on your Kaye Validator AVS or ValProbe RT console and saved on the console and/or shared network drive, can be transferred to the Common Reporting Tool's computer via the sync out functionality for report generation.

The ability to generate reports off-tablet on a separate PC provides added flexibility and efficiency in the use of the Kaye consoles. All files are encrypted, tamper proof, and can only be read by the CRT software, ensuring the same level of data integrity that you receive in our Validation Console, but in the convenience of your office. The interface mirrors the patented Asset Centric Concept included on Kaye Validation Consoles, making for a seamless transition for users. Included in the CRT software is support for Microsoft's Active Directory\*.

Active Directory streamlines user management, ensuring Users are given the correct permissions and access. A free 30-day trial is available for download on the Kaye website. A base license can then be purchased for any number of PCs at your site that would like to view and manage the reporting data.


BENEFITS

- View data and generate reports from a Validator AVS or ValProbe RT system at the convenience of your desk or home office
- Read, manage, and store data on your secured network

- or server, allowing for easy access to files
- Merging of up to three ValProbe RT studies into one single report
- Detailed Pass/Fail Report, including info on passed qualifications
- Report for climatic and thermostatic chambers following AFNOR FD X15-140 regulation for ValProbe RT
- Encrypted and tamper-proof files can only be read through Common Reporting Tool software, ensuring data integrity
- Incorporating Microsoft's Active Directory, ensuring seamless and convenient user management

FD X15-140 Report														
Study Name: Min. load product														
Printed on 21-Mar-2024 at 11:33:36 by JDS														
Temperature Data (°C)														
Sensor Label	SN	Min	Max	Max - Min	Std Dev	Uncertainty	Avg	Uniq	Specification	Avg + Unq	Calc. Spec. Table	Avg - Unq	Calc. Spec. Table	Result
Sensor17	TE18-A	20.96	21.19	0.23	0.06	0.40	21.05	0.42	22.00	21.47	25.00	20.63	20.00	PASS
Sensor12	TI12-A	21.42	21.94	0.52	0.17	0.30	21.70	0.45	22.00	22.15	25.00	21.25	20.00	PASS
Sensor13	U007-A	20.88	21.22	0.34	0.08	0.20	20.04	0.26	22.00	21.20	25.00	20.68	20.00	PASS
Sensor14	U298-B	21.48	21.82	0.34	0.10	0.10	21.63	0.22	22.00	21.85	25.00	21.41	20.00	PASS
Sensor15	M275-B	21.29	21.53	0.25	0.07	0.40	21.39	0.42	22.00	21.81	25.00	20.97	20.00	PASS
Sensor2	BH15-A	20.97	21.95	0.98	0.32	0.40	21.00	0.40	22.00	21.40	25.00	20.80	20.00	PASS
Sensor3	U086-A	20.93	20.99	0.06	0.01	0.40	20.95	0.40	22.00	21.35	25.00	20.55	20.00	PASS
Sensor4	U090-A	21.14	21.24	0.10	0.03	0.40	21.18	0.40	22.00	21.58	25.00	20.78	20.00	PASS
Sensor16	U258-A	21.16	21.26	0.10	0.03	0.20	21.20	0.21	22.00	21.41	25.00	20.99	20.00	PASS
Sensor19	HA19-A	21.21	21.36	0.15	0.04	0.20	21.28	0.22	22.00	21.50	25.00	21.06	20.00	PASS
Sensor5	U052-A	20.83	21.11	0.28	0.06	0.20	20.90	0.23	22.00	21.13	25.00	20.67	20.00	PASS
Sensor6	RP65-A	21.16	21.26	0.12	0.04	0.30	21.22	0.31	22.00	21.53	25.00	20.91	20.00	PASS
Sensor18	PL6-T	20.87	21.38	0.51	0.13	0.30	21.07	0.39	22.00	21.46	25.00	20.68	20.00	PASS
Sensor20	U259-A	21.26	21.46	0.20	0.06	0.30	21.35	0.32	22.00	21.67	25.00	21.03	20.00	PASS
Sensor7	M270-A	21.33	21.45	0.12	0.03	0.30	21.38	0.31	22.00	21.69	25.00	21.07	20.00	PASS
Sensor8	HA18-B	21.35	21.46	0.11	0.03	0.20	21.40	0.21	22.00	21.61	25.00	21.19	20.00	PASS
Sensor9	TV03-A	21.29	21.39	0.10	0.03	0.20	21.33	0.21	22.00	21.54	25.00	21.12	20.00	PASS
Sensor11	U045-A	21.01	21.05	0.04	0.01	0.20	21.03	0.20	22.00	21.23	25.00	20.83	20.00	PASS
Sensor10	U258-B	21.18	21.32	0.14	0.04	0.20	21.24	0.21	22.00	21.45	25.00	21.03	20.00	PASS

AFNOR FD X15-140 Detailed Report



AFNOR FD X15-140 Compliant Graph

Kaye representative contact:

EUROPE, MIDDLE EAST, AFRICA AND ASIA  
Amphenol Advanced Sensors Germany GmbH  
Sinsheimer Strasse 6  
D-75179 Pforzheim  
T: +49 (0) 7231-14 335 0  
F: +49 (0) 7231-14335 29  
Email: kaye@amphenol-sensors.com

USA/AMERICAS  
Amphenol Thermometrics, Inc.  
967 Windfall Road  
St. Marys, PA 15857  
T: +1(814) 834-9140  
F: +1(814) 781-7969  
Email: kaye-us@amphenol-sensors.com

Request a demo:

INDIA  
Amphenol Interconnect India Pvt Ltd.  
Plot no. 6, Survey No.64 | Software Units layout  
MAHAVEER TECHNO PARK  
Hitech City, Madhapur | Hyderabad,  
Telangana – 500081 | T: +91 40 33147100  
Email: kaye-india@amphenol-sensors.com

CHINA  
Amphenol (Changzhou) Connector  
Systems Co., Ltd, Building 10,  
Jintong Industrial Park, No. 8 Xihu Road,  
Wujin High-Tech Development Zone,  
Changzhou, Jiangsu 213164  
T: 0086-519-83055197



**Warranty and disclaimer:** The information mentioned on documents are based on our current tests, knowledge and experience. Because of the effect of possible influences in an application of the product, they do not exempt the user from their own tests, checks and trials. A guarantee of certain properties or a guarantee for the proper suitability of the product for a specific, especially permanent application can not be derived from our data. Liability is therefore excluded to that extent permitted by law. Any proprietary rights of third parties as well as existing laws and regulations must be observed by the recipient of the product on his own responsibility.

© 2024 Amphenol Corporation. All Rights Reserved. Specifications are subject to change without notice. Other company names and product names used in this document are the registered trademarks or trademarks of their respective owners.